

BILI-THERAPY Spot Type

LED Phototherapy Unit



Effective Phototherapy Treatment and Space-saving



Mount Type

Effective phototherapy treatment by blue LEDs

- Promote early hospital discharge by reducing the treatment time in providing high irradiation intensity and effective wavelength.

Space-saving

- The unit can be mounted on the F-Rail without using a stand for space-saving around the incubator or warmer.

Switchable to examination lamp with one-touch of a button

- As an examination lamp, the condition of an infant can be monitored. (White LEDs) The dimmer function that is ideal for developmental care is available.



Examination lamp

Arm Type

Specifications

Atom's code. 60361 **BILI-THERAPY Spot Type** Arm Type
 Atom's code. 60360 **BILI-THERAPY Spot Type** Stand Type
 Atom's code. 60377 **BILI-THERAPY Spot Type** Mount Type

Power requirements	Rating:AC100~240V Frequency:50/60Hz Power consumption:30VA
Classification	Class II, No applied part
Light source	Blue LEDs
Irradiance level	At least30-40 μ W/cm ² /nm (obtained by the BiliBlanket [®] Meter at the irradiation distance of 30cm)
Illumination intensity	Approx.500~2000 lux (Distance from light source:30cm)
Dimensions	Arm Type:100(W) \times 13(D) \times 45(H)cm Stand Type:45(W) \times 71(D) \times 190(H)cm Mount Type (irradiation head):12.5(W) \times 33(D) \times 8(H)cm Mount Type (power source part):14(W) \times 34(D) \times 11(H)cm
Weight	Arm Type:Approx. 2.5kg Stand Type:Approx. 12kg Mount Type (irradiation head):Approx. 0.8kg Mount Type (power source part):Approx. 0.8kg
Accessories	Operation Manual1



Stand Type

Option

Atom's code. 60380 **F-Rail Attachment**
 For BILI-THERAPY Spot Type For Mount Type

Atom's code. 60381 **Pole Mounting Attachment**
 For BILI-THERAPY Spot Type For Mount Type

Manufacturer: ATOM MEDICAL CORPORATION
<http://www.atomed.co.jp>

Exporter: ATOM MEDICAL INTERNATIONAL, INC.
 Iwakata Bldg., 3-18-16, Hongo, Bunkyo-ku, Tokyo, JAPAN
 Phone: + 81 3 3815 2941 Fax: + 81 3 3812 9670/4080
<http://www.atom-ami.co.jp>